



Sepsis Trolley

Sepsis Trolley

Features:

Trolley designed for use when treating Sepsis.

Supports delivery of 'The Sepsis Six' with drawers labelled to allow completion of the BUFALO bundle within 1 hour of sepsis recognition

- Administer Oxygen
- Take Blood Cultures
- Give IV Antibiotics
- Give IV Fluids
- Check Venous Blood Gas Lactates
- Measure Urine Output

Drawer locking methods:

- Battery-operated push-button lock
- High security bolt lock – supplied with two keys
- Pull tight seal and full-length bar

5 labelled drawers (IV Fluids and IV Antibiotics in same drawer)

- Oxygen
- Blood Cultures
- IV Fluids & Antibiotics
- Lactates
- Urine Output

Supplied accessories

- Push handle
- 2 x divider sets (FOC)
- Side pull out extension table

Available accessories

- Oxygen tank holder
- Glove box holder
- Sharps bin holder
- IV pole
- Other accessories can be provided upon request.



Oxygen tank holday



Glove box holder



IV pole



Overall dimensions:

Agile Medical sepsis trolley is available in two sizes, **narrow** and **standard**.

Narrow dimensions (mm): Width 520, Depth 540, Height 1046

Standard dimensions (mm): Width 690, Depth 540, Height 1110

Standard colours:

Trolley sides – orange | **Drawer front** – white | **Work surface** – white.

Options

Sharps bin holder

SKU	Description	Storage Capacity
SBH100	Sharps Bin Holder with Strap	<5L
SBH200	Sharps Bin Holder	<7L
SBH300	Sharps Bin Holder	<14L

- All sharps bin holders secured to trolley side rail
- With SBH100, sharps bin secured using Velcro strap
- With SBH200 and SBH300, sharps bin stored in wire basket



Universal sharps bin mountable bracket

SKU	Description	Storage Capacity
SSMB01	Sharps Bin Mountable Bracket	<32L

- Mounted to trolley side rail
- Compatible with following sharps bins: Daniels Healthcare S14, S22 and S32



Security Device

SKU	Description	Length
DRS200	Security Device	1800mm

- Allows the trolley to be secured to the 'fabric' of the building
- Additional keys provided
- The fixed end is secured to the wall, the free end is locked to the steel corners of the trolley when not in use
- Recommended wall fitting height 1000mm



Important Note - Security Clamps

It is a requirement that when the trolley is not in use it should be secured to the "fabric of the building". Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area and broken into, but also encourages good nursing practice, giving the trolley a "permanent home".

"Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB- March 2005"



Head Office

Agile Medical
Estate Road 2
South Humberside Industrial Estate
Grimsby
DN31 2TG
UK

London Office

Unit 1 Kimberly Business Park
Blackness Lane
Keston
BR2 6HL
UK

Phone: +44 (0)1472 266 190

Email: sales@agilemedical.co.uk

Web: agilemedical.co.uk

Agile
M E D I C A L